




RENTAL SOLUTIONS

# Rental product portfolio

# NEWER RENTAL FLEET

**95%** is less than 5 years of age

Arjo is consistently investing in our fleet and processes to ensure that you always get the most up-to-date equipment possible.



## Medical beds



### Total Lift™ Bed

Total Lift Bed can progressively move patients to full standing and weight bearing to aid rehabilitation early in the acute phase of critical illness.

- Weight limit: 425 lbs



### Enterprise® 9000 hospital bed

Universal hospital bed offering solutions for the daily needs of caregivers and patients. Suitable for low to high acuity.

- Weight limit: 407 lbs



### Citadel™ Plus bariatric care system

Bariatric bed designed to deliver value and dignified care for patients of size with flexible and safe solutions.

- Weight limit: 1,000 lbs



### Citadel™ C200 integrated system

Flexible, integrated therapeutic mattress and patient care system for high-risk patients. Suitable for medurg, critical care and ICU.

- Weight limit: 500 lbs

# Therapeutic support surfaces

We offer a comprehensive portfolio of therapeutic support surfaces designed to manage patient outcomes by helping to prevent complications associated with immobility. From reactive to active surfaces, our portfolio includes hybrid, low air-loss, alternating pressure surfaces and microclimate management to help prevent the various stages of pressure injuries.



## AtmosAir® Velaris adaptable alternating system with or without pump

- Adaptable alternating pressure with patented air redistribution module to maintain a predefined level of pressure in reactive mode, and AltoVac vacuum technology to provide full offloading in active mode. Can instantly step up/down therapy with addition of pump
- Weight limit: 550 lbs



## AtmosAir® Plus bariatric hybrid

- Bariatric pressure redistribution surface
- Weight limit: 1,000 lbs



## First Step® Cirrus low air-loss

- Low air-loss technology with high performance medical fabric and non-slip edging
- Weight limit: 300 lbs



## MaxxAir® ETS bariatric low air-loss

- Bariatric low air-loss therapy with full body turn assist up to 30 degrees, continuous lateral rotation and excellent MVTR
- Weight limit: 1,000 lbs



## Auralis® automatic alternating system

- Active therapy (alternating) & reactive therapy (constantly low pressure modes) with intelligent Self-Set Technology that automatically senses patient weight to adjust pressure settings accordingly
- Weight limit: 550 lbs



## Auralis® Plus bariatric automatic alternating system

- Bariatric active therapy (alternating) & reactive therapy (constantly low pressure modes) with intelligent Self-Set Technology that automatically senses patient weight to adjust pressure settings accordingly
- Weight limit: 1,000 lbs



## Nimbus® Pro alternating pressure surface

- Alternating pressure surface that uses vent valve technology to selectively isolate and deflate any one of the 19 individual cells under the patient's body for customized off-loading
- Weight limit: 550 lbs



## Skin IQ® Microclimate Management

- Waterproof, vapor-permeable top layer that offers a bacterial and viral barrier<sup>1</sup>
- Negative airflow Technology (NAT) continually draws excess moisture away from the skin/mattress interface

### Reference

1. Arjo Data on File. Test Report Viral Barrier Property of Top Sheet of Skin IQ Coverlet — Test Report 100018376 pg. 3 & Test Report 100019009 pg. 1-7.

# Bariatric accessories solutions



## Wheelchair

Sturdy oversized steel frame to safely transport patients up to 700 lbs (317.5 kg)



## Walker

Oversized, anodized aluminum, double-tubed frame — strong enough to safely support patient weighing up to 850 lbs (385.5 kg)



## Commode/shower chair

Oversized, anodized aluminum, double-tubed frame — strong enough to safely support patient weighing up to 1,000 lbs (453.5 kg)



# Critical care solutions

Arjo provides a range of solutions that allow clinicians to safely and efficiently reposition critically ill, immobile patients. We offer critical care solutions that help our customers **minimize risk and maximize efficiency**.



## RotoProne®

Comprehensive system that provides an efficient and effective way to automatically prone patients.

- Weight limit: 88-350 lbs
- Automated prone therapy



## TriaDyne Proventa®

A comprehensive system combining pulmonary and kinetic therapies clinically proven to improve pulmonary outcomes and maintain skin integrity associated with critically ill, immobilized patients.

- Weight limit: 350 lbs
- Percussion therapy



## RotoRest® Delta

The RotoRest Delta provides advanced kinetic therapy to mobilize the most immobile patients who are critically ill without compromising their condition.

- Weight limit: 300 lbs
- Kinetic therapy — continuous lateral rotation therapy (CLRT) up to 62°



## FluidAir® Elite

Provides full body air-fluidized therapy for advanced pressure redistribution and microclimate management, helping to heal skin and relieve pain for patients with pressure injuries, burns, surgical flaps and grafts, and other skin wounds.

- Weight limit: 250 lbs
- Air-fluidized therapy system



# Clinical, therapeutic seating



## Phoenix™ chair

- Weight limit: 364 lbs
- Suitable for highly dependent patients with complex needs



## Milano™ chair

- Weight limit: 364 lbs
- Designed for accessibility, it can easily be adjusted to fit different patients and residents



## Atlanta™ chair

- Weight limit: 364 lbs
- A robust and durable chair, it provides a safe and comfortable seat for those with vigorous involuntary movements, while reducing the risk of falling from the chair.



## Sorrento™ chair

- Weight limit: 364 lbs
- Designed to maximize comfort, it provides postural support and pressure redistribution. The tilt-in-space chair is clinically designed for the patient and caregiver in mind



## Bariatric Sorrento™ chair

- A tilt-in-space bariatric chair on wheels and can support patients up to 650 lbs
- Designed to reduce the risk for injury. Fully motorized to limit patient handling, accommodates unique body shapes, limits sliding, provides postural support and aids in repositioning.



# All of our rental equipment is available for you to order online through our **Arjo Express Portal**.

Taking care of your patients with a pressure injury and critical needs is easier than ever before. Set up online management of therapeutic surfaces and medical beds for your facility.

## Arjo Express offers many benefits, including:

- Secure and confidential site access
- Tracking of Arjo therapeutic surface systems throughout one or multiple facilities
- Valuable rental information including: rental start date, rental end date, patient identifier, unit serial number, and unit location
- Rental reporting, including daily, monthly, quarterly and annually
- Ease of unit-to-unit transfers

Phone  
**+1-800-323-1245**  
24hr Rental hotline  
**+1-800-343-0974**  
or visit  
**Arjo.us/Rental**

## Product resources



**Enterprise 9000**



**Citadel**



**Citadel Plus**



**Total Lift Bed**



**AtmosAir Velaris**



**Auralis**



**First Step Cirrus**



**Skin IQ**



**MaxxAir ETS**



**Sorrento**



**Bariatric Sorrento**



**Auralis Plus**

At Arjo, we believe that empowering movement within healthcare environments is essential to quality care. Our products and solutions are designed to promote a safe and dignified experience through patient handling, medical beds, personal hygiene, disinfection, diagnostics, and the prevention of pressure injuries and venous thromboembolism. With over 6500 people worldwide and 65 years caring for patients and healthcare professionals, we are committed to driving healthier outcomes for people facing mobility challenges.

Arjo Inc., 2349 West Lake Street · Addison, Illinois 60101 · 1-800-323-1245 · [www.arjo.us](http://www.arjo.us)

